



Cosmopor® E sterile adhesive, island wound dressing

For sterile dressing of post operative wounds and/or minor injuries



Viscose Core

- Provides high absorbency and cushioning

Polyethylene net contact layer

- Prevents adherence of the dressing to the wound

Hypoallergenic fibres, permeable to air and water

- Prevents irritation of surrounding skin
- Enables the natural functions of the skin to continue



PAUL HARTMANN LTD
 Heywood Distribution Park
 Pilsworth Road
 Heywood
 Lancashire
 OL10 2TT
 Tel: 01706 363200
 Fax: 01706 363201

Web: www.hartmann.co.uk
 E-mail: info@uk.hartmann.info

Cosmopor E sterile adhesive, island wound dressing

For sterile dressing of post operative wounds and/or minor injuries

Cosmopor E is a sterile, self-adhesive wound dressing made of soft non woven Polyester, with an absorbent viscose pad covered with a non adherent porous wound contact layer. This is coated with skin friendly synthetic adhesive.

Cosmopor E carries the CE mark according to EU directive 93/42/EEC for medical devices. The product is classified as a class I sterile medical device. A conformity assessment has been performed for Cosmopor E and it has been shown to be in compliance with all applicable requirements of the above mentioned directive.

Indication

Cosmopor E is indicated for post operative wound management and sterile dressing of minor injuries e.g. in First Aid.



Product	Size	Pack	Hartmann Code	NHS Code	PIP Code
Cosmopor® E	7.2 x 5cm	50 pcs	900 990	EIJ038	273 1222
	10 x 8cm	25 pcs	900 993	EIJ039	273 1255
	15 x 8cm	25 pcs	900 994	EIJ040	273 1263
	20 x 10cm	25 pcs	900 996	EIJ041	273 1289
	25 x 10cm	25 pcs	900 997	EIJ042	273 1297
	35 x 10cm	25 pcs	900 998	EIJ043	273 1305