



# Hospicrepe® Cotton Crepe Bandage

For light support and fixation



### 100% Cotton with 120% stretch capacity

- Ideal for all general bandage requirements
- Two varieties available: 233 and 239

### Woven edges

- Reduces bandage fray



**PAUL HARTMANN LTD**  
 Heywood Distribution Park  
 Pilsworth Road  
 Heywood  
 Lancashire  
 OL10 2TT  
 Tel: 01706 363200  
 Fax: 01706 363201

**Web: [www.hartmann.co.uk](http://www.hartmann.co.uk)**  
**E-mail: [info@uk.hartmann.info](mailto:info@uk.hartmann.info)**

# Hospicrepe® Cotton Crepe Bandage

## For light support and fixation



Hospicrepe® is a cotton crepe bandage ideal for all general bandage requirements including support, pressure and fixation. The crepe character is achieved with highly twisted 100% cotton yarn in the weft and non-elastic yarn in the warp.

Two varieties available: 233 BP Standard and 239 (hospital quality) which has a thinner yarn and less number of threads for greater economy.

### Indication

Hospicrepe® is indicated for support and relief bandaging of distortion and dislocations and for fixation.



Product	Variant	Size	Pack	Hartmann Code	NHS Code
Hospicrepe® cotton crepe bandage	233	5cm x 4.5m	12 pcs	915 561	ECA 811
	233	7.5cm x 4.5m	12 pcs	915 562	ECA 030
	233	10cm x 4.5m	12 pcs	915 563	ECA 009
	233	15cm x 4.5m	12 pcs	915 564	ECA 059
	239	5cm x 4.5m	12 pcs	915 557	ECA 088
	239	7.5cm x 4.5m	12 pcs	915 558	ECA 089
	239	10cm x 4.5m	12 pcs	915 559	ECA 090
	239	15cm x 4.5m	12 pcs	915 560	ECA 091