



Parklands Distribution Centre  
Pilsworth Road  
Heywood  
Lancashire OL10 2TT  
Telephone: 01706 363200

## Hydrosorb Fluid Handling Capacity

<b>Time Interval</b>	<b>Absorbency g/10cm<sup>2</sup></b>	<b>MVTR g/10cm<sup>2</sup></b>	<b>Fluid Handling Capacity g/10cm<sup>2</sup></b>
24 hours	1.28	0.31	1.59
48 hours	1.17	0.51	1.69
72 hours	1.14	0.77	1.91
96 hours	1.07	1.08	2.15