



helps healing.

# Omnistrip® Adhesive Wound Closure Strips

For atraumatic closure of minor wounds and surgical incisions



#### Hypoallergenic, pressure sensitive adhesive

- Ensures reliable and secure wound closure
- Prevents irritation of surrounding skin
- Can be removed without pain or residue

#### Microporous, neutral coloured non woven fabric

- Avoids maceration and allows the skins natural functions to continue



HARTMANN



helps healing.

**PAUL HARTMANN LTD**

Heywood Distribution Park  
Pilsworth Road  
Heywood  
Lancashire  
OL10 2TT  
Tel: 01706 363200  
Fax: 01706 363201

Web: [www.hartmann.co.uk](http://www.hartmann.co.uk)

E-mail: [info@uk.hartmann.info](mailto:info@uk.hartmann.info)

# Omnistrip® Adhesive Wound Closure Strips

## For atraumatic closure of minor wounds and surgical incisions

Omnistrip® is a wound closure strip made of neutral coloured, radio translucent non woven fabric. The hypoallergenic polyacrylate adhesive means Omnistrip® has low irritant properties and can be removed without pain or residue.

Omnistrip® is comfortable to wear with rounded edges, temperature resistant and is well tolerated by the skin due to extremely high permeability to air and water vapour.

Using Omnistrip® instead of sutures or staples will reduce scarring.

### Indication

For atraumatic closure of minor wounds and surgical incisions; to protect wound margins against tension; to stabilise fresh scars and to provide wound support following early suture or staple removal.

Ensure that skin is clean and dry before application.



Product	Size	Strips per envelope	Envelopes per pack	Hartmann Code	PIP Code
<b>Omnistrip® Wound closure strip</b>	3 x 76mm	5 pcs	50	540681	—
	6 x 38mm	6 pcs	50	540682	—
	6 x 76mm	3 pcs	50	540683	353 1944
	6 x 101mm	10 pcs	50	540684	—
	12 x 101mm	6 pcs	50	540685	—
	25 x 127mm	4 pcs	50	540686	—