

Peha-haft® Latex Free Cohesive Conforming Bandage

For safe and effective dressing retention



Latex Free, first of it's kind!

- Specially designed for patients prone to allergies and persons allergic to latex
- Optimum in functionality and user friendliness
- Odourless

Breathable fabric with double adhesive effect

- High proportion of natural fibres
- Permeable to air
- Very kind to skin due to the absence of acrylate and chlorine



PAUL HARTMANN LTD
 Heywood Distribution Park
 Pilsworth Road
 Heywood
 Lancashire
 OL10 2TT
 Tel: 01706 363200
 Fax: 01706 363201

Web: www.hartmann.co.uk
E-mail: info@uk.hartmann.info

Peha-haft® Latex Free Cohesive Conforming Bandage

For safe and effective dressing retention

Peha-haft and Peha-haft® latex free are cohesive, elastic conforming bandages specifically designed to facilitate all kinds of dressing retention. Breathable fabric with a special double adhesive effect means that the bandage sticks to itself, but not to skin.

85% extensibility with a high absorption capacity due to a high number of natural fibres, the bandage consists of 40% Cotton, 30% Viscose and 30% Polyamide

Benefits:

- Provide safe dressing retention even on joints and conical or round body parts.
- No additional fixation of the bandage needed.
- Cost saving due to low consumption of material.


There is only one thing you can do against latex allergy: Avoid Latex

In Europe about 30% of the population suffer from allergies. The immediate allergy to natural latex rubber is widely spread and presents a considerable problem in health care.

Peha-haft® latex free offers perfect protection for staff and for patients.

Indication

For fixation bandaging, especially on joints and conical or round parts of the body; retention of wound dressings and of orthopaedic padding material, cannulae, etc and for use as support bandage providing light compression.



Source: European Centre for Allergy Research Foundation – ECARF

Product	Size	Pack	Hartmann Code	NHS Code	PIP Code
Peha-haft® Latex Free cohesive, conforming bandage	4cm x 4m	single	932 441		347 3956
	6cm x 4m	single	932 442	EBA 054	347 3964
	8cm x 4m	single	932 443	EBA 055	347 3972
	10cm x 4m	single	932 444	EBA 056	347 3980
	15cm x 4m	single	932 445		347 3998
Peha-haft®	4cm x 4m	single	932 421		
	6cm x 4m	single	932 423		
	8cm x 4m	single	932 424		
	10cm x 4m	single	932 425		
	10cm x 20m	single	932 435		
	12cm x 4m	single	932 426		