



# Sterifix<sup>®</sup> sterile absorbent dressing with adhesive strips

For dressing of simple wounds with slight exudation



#### Cotton wool core

- Provides skin-friendly absorbency

#### Microporous, polyethylene contact layer

- Prevents adherence to wound

#### PAUL HARTMANN LTD

Heywood Distribution Park  
Pilsworth Road  
Heywood  
Lancashire  
OL10 2TT  
Tel: 01706 363200  
Fax: 01706 363201

Web: [www.hartmann.co.uk](http://www.hartmann.co.uk)

E-mail: [info@uk.hartmann.info](mailto:info@uk.hartmann.info)

Product	Size	Pack	Hartmann Code	NHS Code	PIP Code
Sterifix <sup>®</sup>	5 x 7cm	10 pcs	241 950	–	266 7319
	7 x 10cm	10 pcs	241 951	–	266 7327
	10 x 14cm	10 pcs	241 952	–	266 7335