



# Hospilite<sup>®</sup> Permanently Elastic Crepe Bandage

For light support and fixation

**TYPE 2**  
LIGHT SUPPORT



**96% Cotton, 3% Polyamide, 1% Elastane**

- Ideal for all general bandage requirements
- Easy to apply

**Woven edges**

- Reduces bandage fray



**PAUL HARTMANN LTD**  
 Heywood Distribution Park  
 Pilsworth Road  
 Heywood  
 Lancashire  
 OL10 2TT  
 Tel: 01706 363200  
 Fax: 01706 363201

**Web: [www.hartmann.co.uk](http://www.hartmann.co.uk)**  
**E-mail: [info@uk.hartmann.info](mailto:info@uk.hartmann.info)**

# Hospilite® Permanently Elastic Crepe Bandage

## For light support and fixation



**TYPE 2**  
**LIGHT SUPPORT**



Hospilite® is a cost effective permanently elastic crepe bandage, ideal for all general bandage requirements including support, pressure and fixation.

### Indication

Hospilite® is indicated for support and relief bandaging of distortion and dislocations and for fixation.

A Type 2 Light Support Bandage is described in BS7505:1995 as a 'Bandage used to restrict movement and to provide intermittent sub-bandage pressure'.

Product	Size	Pack	Hartmann Code	NHS Code	PIP Code
Hospilite® permanently elastic crepe bandage	5cm x 4.5m	12 pcs	935151	ECA193	347 0192
	7.5cm x 4.5m	12 pcs	935152	ECA194	347 0200
	10cm x 4.5m	12 pcs	935153	ECA195	347 0218
	15cm x 4.5m	12 pcs	935154	ECA196	347 0226